61-5005 / 61-5007

H/S Catheter Set (English)

H/S Katheter-Set (Deutsch / German)

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Directions for Use

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Mode d'emploi

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61-5005 and 61-5007 • H/S Catheter Set for Hysterosonography and Hysterosalpingography

Directions for Use (English)

DEVICE DESCRIPTION

The H/S Catheter Set consists of a latex-free balloon catheter, insertion sheath and a 1.5 cc (ml) syringe (5F catheter) or 3.0 cc (ml) syringe (7F catheter). The catheter can be used with aqueous based contrast media. (See *Figure 1.*)


Figure 1

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE/INDICATIONS

For administering contrast media during Hysterosonography or Hysterosalpingography procedures to detect uterine pathology such as polyps, fibroids, adhesions or endometrial thickening, or patency of fallopian tubes.

CONTRAINDICATIONS

Suspected infection, suspected pregnancy, profuse bleeding or sexually transmitted disease.

PRECAUTIONS

Do not exceed the recommended balloon inflation volume of 1.5 cc (ml) (5F catheter) or 3.0 cc (ml) (7F catheter) or the balloon may burst.

The use of OIL-BASED contrast media such as ethyl esters may interact with the balloon of the catheter, causing rupture. The use of oil-based contrast media is not recommended.

ADVERSE EVENTS

Some patients may have a hypersensitivity to contrast media.

DIRECTIONS FOR USE

Catheter Preparation

1. Grasp the translucent insertion sheath connected to the trapezoidal hub.
2. Remove and discard the crimped yellow protective cover, exposing the balloon catheter tip.
3. Test the balloon integrity by inflating with air, saline or water using the syringe in the set. If performing Hysterosonography remove as much air as possible from the balloon. Deflate *completely* by pulling back on the syringe plunger and closing the stopcock.
4. Attach a contrast media-filled syringe (not supplied) to the blue luer hub and fill the catheter with contrast media to expel air. Advance the insertion sheath so that the distal end of the catheter protrudes slightly from the distal end of the sheath. (See *Figure 2.*)

Figure 2

Catheter Placement

5. Visualize the external cervical os with the aid of a speculum and advance the sheath and catheter so that the tip of the catheter enters the cervical canal.
6. Hold the sheath stationary and advance the catheter into the cervical canal and into the uterine cavity.
7. Open the stopcock and **slowly** inflate the balloon with up to 1.5 cc (ml) (5F catheter) or 3.0 cc (ml) (7F catheter) of air, saline or water. Turn the stopcock off allowing the balloon to remain inflated.
8. Withdraw the catheter until it occludes the internal os of the cervix.

Injection of Contrast Media

9. Inject the contrast media, close the clamp and complete the study in a routine manner.
10. Occasionally, when access to the uterine cavity is difficult the balloon must be inflated within the endo-cervical canal.

Catheter Removal

11. Open the stopcock and deflate the balloon by pulling back on the syringe plunger. Withdraw the catheter.

EXPLANATION OF SYMBOLS

Reorder Number

Batch Code

Use By

Irradiation Sterilization

Do Not Reuse

Latex Free

ATTENTION: See instructions for use.

Product conforms to the Medical Device Directive 93/42/EEC

Authorized Representative in the European Community.

95 Corporate Drive • Trumbull, CT 06611 USA
Phone: (800) 243-2974 • Fax: (800) 262-0105

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© 2011 CooperSurgical, Inc.
Leisegang Feinmechanik GmbH