ASPIRATORE
ASPIRATOR
ASPIRATEUR
ABSAUGER
ASPIRADOR

NEW **ASKIR** 30

30751/614 - Rev. 7 del 17.09.2012

Italian
Medical
Touch

MANUALE D'USO

INSTRUCTION MANUAL

MODE D'EMPLOI

HANDBUCH

MANUAL DE ISTRUCCIONES

CE 0123

CA-MI srl - Via Ugo La Malfa 13 - 43010 Pilastro (Pr) - Italy
Tel. +39 0521 637133 +39 0521 631138 - Fax +39 0521 639041
export@ca-mi.it - vendite@ca-mi.it - www.ca-mi.it

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NEW ASKIR 30 Surgical aspirator is a portable unit, working with 230V ~ / 50 Hz network electricity, designed for the aspiration of bodily fluids in adult and children. It's particularly suitable for nasal, oral or tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another. Easily portable equipment designed for continuous use.

Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European

safety standard, the product is supplied with a complete polycarbonate autoclavable jar with overflow valve and it is equipped with aspiration regulator and vacuum indicator located on the front panel.

GENERAL WARNING

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

ONLY HIGHLY QUALIFIED STAFF USE RESERVED

THE INSTRUMENT MUST NOT BE DISASSEMBLED. FOR A TECHNICAL SERVICE ALWAYS CONTACT CA-MI

IMPORTANT SAFETY RULES

1. Check the condition of the unit before each use. The surface of the unit should be carefully inspected for visual damage. Check the mains cable and **do not connect to power** if damage is apparent;
2. Before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Respect the safety regulations indicated for electrical appliances and particularly:
 - Use original components and accessories provided by the manufacturer CA-MI to guarantee the highest efficiency and safety of the device;
 - The device can be used only with the bacteriological filter;
 - Never immerse the appliance into water;
 - Place instrument on stable and flat surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Don't use in the presence of inflammable substances such as anaesthetic, oxygen or nitrous oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - Keep off the reach of children or not capable people without supervision;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
 - Don't use the device thoracic drainage.
5. For repairs, exclusively contact CA-MI technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer will not be responsible for damage due to improper use or connection to an electrical system not complying with current regulations;
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
8. Instrument and accessories discharging must be done following current law regulations in every country of use.
9. None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact CA-MI technical assistance
10. Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same.

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

In respect of art. 13 Decreto Legislativo 25 Luglio 2005, n.151 "Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal"

The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health. In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.

TECHNICAL CHARACTERISTICS

TYPOLOGY (MDD 93/42/EEC)	Dispositivo Medico Classe IIa
MODEL	NEW ASKIR 30
CLASSIFICATION UNI EN ISO 10079-1	HIGH VACUUM / HIGH FLOW
MAIN VOLTAGE	230V~ / 50Hz
POWER CONSUMPTION	107VA
FUSE	F 1 x 1.6A L 250V
MAXIMUM SUCTION PRESSURE (without jar)	-80kPa (-0.80 Bar)
MINIMUM SUCTION PRESSURE (without jar)	Less - 40kPa (-0.40 bar)
MAXIMUM SUCTION FLOW (without jar)	40 l/min
WEIGHT	3.6 Kg
SIZE	350 x 210 x 180 mm
DUTY CYCLE (to 35°C and 110% operating voltage)	Non-Stop Operated
SICILICONE TUBE SIZE	Ø 6 x 10 mm
ACCURACY OF VACUUM INDICATOR	± 5%
WORKING CONDITION	Room temperature: 5 ÷ 35°C Room humidity percentage: 30 ÷ 75% RH Atmospheric pressure: 700 ÷ 1060 hPa
CONSERVATION CONDITION AND TRASPORT	Room temperature: -40 ÷ 70°C Room humidity percentage: 10 ÷ 100% RH

SYMBOLS

	Class II isolation equipment
CE 0123	CE marking in conformity with EC directive 93/42/EEC and subsequent changes
	General warnings and/or specifications
	Consult the instruction manual
	Manufacturer: CA-MI S.r.l. Via Ugo La Malfa nr.13 – 43010 Pilastro (PR) Italia
	Keep in a cool, dry place
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
Hz	Mains Frequency
I	ON
0	OFF
LOT	Lot Number
SN	Serial Number
REF	Identification device
IP21	Degree of protection an electrical device provides in the case of accidental or intentional contact with the human body or with objects, and protection in the case of contact with water.

Please note technical specifications may vary upon the manufacturer's discretion!

This section contains information regarding the conformity of the compliance with the IEC 60601-1-2 Standard.

The NEW ASKIR 30 surgical aspirator is an electro-medical device that requires particular precautions regarding electro-magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied. Mobile and portable RF communication appliances (mobile phones, transceivers, etc..) can affect the medical system.

The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system.

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The surgical aspirator NEW ASKIR 30 is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator NEW ASKIR 30 should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator NEW ASKIR 30 only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances. The surgical aspirator NEW ASKIR 30 can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Irradiated / Conducted emissions CISPR11	Class [B]	
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Immunity Emissions			
The surgical aspirator NEW ASKIR 30 is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator NEW ASKIR 30 should assure that it's used in such an environment.			
Immunity Test	Level indicated by the IEC 60601-1-2	Compliance Level	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	$\pm 6\text{kV}$ on contact $\pm 8\text{kV}$ in air	The device doesn't change its state	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	$\pm 2\text{kV}$ power supply lines $\pm 1\text{kV}$ for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Surge IEC/EN 61000-4-5	$\pm 1\text{kV}$ differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	$5\%U_T$ ($>95\%$ dip U_T) for 0.5 cycle $40\%U_T$ ($>60\%$ dip U_T) for 5 cycle $70\%U_T$ ($>30\%$ dip U_T) for 25 cycle $<5\%U_T$ ($>95\%$ dip U_T) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator NEW ASKIR 30 request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Note U_T is the value of the power supply voltage			

Guidance and manufacturer's declaration – Immunity Emissions			
The surgical aspirator NEW ASKIR 30 is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator NEW ASKIR 30 should assure that it's used in such an environment.			
Immunity Test	Level indicated by the IEC 60601-1-2	Compliance level	Electromagnetic environments - guidance
Conducted Immunity IEC / EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	The portable and mobile RF communication devices, including cables, must not be used closer to the NEW ASKIR 30 device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance
Radiated Immunity IEC / EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	$E_1 = 3 \text{ V / m}$	$d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P}$ from 80 MHz to 800MHz $d = [23 / E_1] \sqrt{P}$ from 800 MHz to 2.5 GHz Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site ^{a)} , could be lower than the level of conformity of each frequency interval ^{b)} . It is possible to check for interference in proximity to devices identified by the following symbol: 
Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by reflection from buildings, objects and people. a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary. b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.			

Recommended separation distance between portable and mobile radio-communication devices and the monitor			
The NEW ASKIR 30 surgical aspirator is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the NEW ASKIR 30 device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the NEW ASKIR 30 device, as recommended below, in relation to the radio-communication maximum output power.			
Maximum nominal output power of the Transmitter W	Separation distance from the frequency transmitter (m)		
	150 kHz to 80 MHz $d = [3.5 / V_1] \sqrt{P}$	80 MHz to 800 MHz $d = [12 / E_1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23 / E_1] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.			
Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied			
Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.			

ACCESSORIES SUPPLIED

DESCRIPTION	CODE
COMPLETE ASPIRATION JAR 1000cc	RE 210001
CONICAL FITTING	25723
TUBES SET 6 mm x 10 mm	RE 210410
ASPIRATION PROBE CH20	51100
ANTIBACTERIAL FILTER AND HYDROPHOBIC (+ n°1 spare parts)	SP 0046

Available under request with different versions with complete jar 2000cc.

Antibacterial Filter: The filter is produced with (PTFE) hydrophobic material which prevents fluids entering the pneumatic circuit. The filter is for a single patient use which will protect patients and machines from cross contamination. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately. In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted.

Suction catheter: Single-use device to be used on a single patient. Do not wash or re-sterilize after use. Reuse may cause cross-infections. Don't use after lapse of the sell-by date

WARNING: Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

CLEANING OF ACCESSORIES

To clean the plastic housing of the device wear disposable latex gloves and clean with denaturated alcohol or hypochlorite solutions. Washing and / or cleaning the autoclavable jar as to be carried out as follows:

- Wear protection gloves and apron (glasses and face mask if necessary) to avoid contact with contaminating substances;
- Disconnect the jar from the device
- Disconnect all tubes from the jar and the protection filter
- Empty and dispose of the content and of the suction catheter according to the laws in force in your country;
- Separate all parts of the cover (overflow valve, o-ring);

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly.

Then soak in warm water (temperature shall not exceed 60°C).

Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations.

Rinse in running warm water and dry all parts with a soft cloth (non-abrasive).

The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure – 15 min) making sure that the jar is positioned upsidedown.

Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1). Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended.

After sterilization and cooling at environment temperature of the parts make sure that these are not damaged.

Assemble the jar as follows:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sure that cover seals perfectly to avoid vacuum leakages or liquid exit

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C (1 bar relative pressure – 15 min)

The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C. (1 bar relative pressure – 15 min)

The device is ready for a new employment now.

DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER

MAINTENANCE

The **NEW ASKIR 30** suction equipment does not need maintenance or lubrication.

It is, however, necessary to inspect the unit before each use. Unpack the instrument and **always check** the plastic parts for any damage that may have occurred during prior use. Connect cable to electrical network and turn switch on.

Close the aspiration outlet with your finger and with suction regulator at maximum check that the vacuum position reaches at least -80 kPa (-0.80 bar) maximum. Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down to -40 kPa (-0.40 bar). Check that no loud noises are present.

A fuse (**F 1x1.6A L 250V**) is located in the plug to protect the instrument. If required, replace with another or the correct type.

Fault type	Cause	Solution
1. The suction unit doesn't work	Cable is damaged	Replace the cable
	External power source failure	Check the external power source
2. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
3. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
4. The Vacuum power on the patient side is either very low or absent	a) Vacuum regulator set to minimum b) Protection filter blocked or damaged c) Connection tubes blocked, kinked or disconnected d) Shut-off valve blocked or damaged e) Pump motor damaged	a) Turn the vacuum regulator clockwise and check the value of the vacuum on the gauge b) Replace the filter c) Replace or reconnect the tubes, check the jar connections d) Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position e) Refer to authorised service personnel
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into its place
6. The float doesn't close	The float is covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
Faults 1 - 2 - 3 - 4 - 5 - 6 - 7	None of the remedies has achieved the desired results	Contact the seller or CA-MI After-sales Assistance Service

If the overfill security system fails, continued aspiration will be stopped by the bacteriological filter.

Should this back-up security also fail, there is a risk of penetration of liquid into the device. Don't attempt to proceed with aspiration under these circumstances.

Return the device to CA-MI for service.

CA-MI Srl will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT CA-MI TECHNICAL SERVICE. THE MANUFACTURER DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

CLEANING DEVICE

Use a soft dry cloth with not – abrasive and not – solvent detergents.


BEFORE STRATING CLEANING OPERATIONS, MAKE SURE THAT INTERNAL PARTS OF DEVICE DON'T COME INTO CONTACT WITH LIQUIDS

INSTRUCTIONS

- Place the unit on a flat, horizontal surface
- Connect one end of the short silicon tube, with antibacterial filter, to the suction connector on the lid of the jar.
- The other tube already connected to the filter has to be connected to the "VACUUM" jar outlet, where has been fixed the red float (security float). When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device.

WARNING: Ensure that the FLUID SIDE or IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM". A wrong connection causes immediate destruction in case of contact with sucked liquids.

Filter assembling

- Connect the long silicone tube to the "PATIENT" jar outlet
- Connect the other end of the long silicon tube to the probe plastic connector, then connect the suction probe to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch on position I to start suction.
- Unscrew the lid of the jar and fill the jar 1/3 full of ordinary water (this assists the unit to reach peak vacuum performance and makes clean-up easier) then re-screw the lid on the jar correctly.
- During operation the jar has to be in vertical position to avoid overflow valve to cut off aspiration. Should this happen, switch off the device and disconnect the tube from the jar cover (from "VACUUM" outlet).
- Once finished push switch on 0 position and unplug.
- Remove the accessories and clean.

NEVER USE THE DEVICE WITHOUT JAR AND / OR PROTECTION FILTER

MAKE SURE THAT CHILDREN AND/OR MENTALLY ILL PEOPLE DO NOT USE THE DEVICE WITHOUT ADULT SURVEILLANCE

RULES FOR RETURNING AND REPAIRING

COMPLYING WITH THE NEW EUROPEAN RULES, CA-MI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE.

THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.

CA-MI warrants its products for **24 months** after purchasing date.

In front of this warranty, CA-MI will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective.

The product must be accompanied by a description of the defect.

The warranty, with exclusion of responsibility for direct and indirect damages, it is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of CA-MI for damages caused by the transport or dismay from the vector.

Every returned instrument will be hygienically checked before repairing. If CA-MI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter.

CA-MI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, CA-MI will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

CA-MI is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting.

We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that CA-MI can establish whether it falls into the category of the faults covered by the guarantee.

CA-MI Srl cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.

Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC (and subsequent changes) and its normatives.

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